

Fiberoptic bronchoscopy



Fiberoptic bronchoscopy


- Diagnostic bronchoscopy
 - Evaluate airway anatomy
 - Identify specific etiology of pneumonia
- Therapeutic bronchoscopy
 - Restoration of airway patency
 - BAL in patient with alveolar filling disorder
 - Special procedures: transbronchial biopsy, bronchoscopic intubation, administer drug

Fiberoptic Bronchoscope

- 3 standard pediatric size
 - 3.5 mm OD with 1.2 mm suction
 - 2.8 mm OD with 1.2 mm suction
 - 2.2 mm OD without suction (ultrathin)



Light source & monitoring



Setting- bedside

